



BFT- 15 Portable Medical Refrigerator Unit with Cart

The BFT-15 portable medical BioFridge with cart was designed for use within hospitals, laboratories and medical facilities.

The internal power supply will run all units for up to (18) eighteen hours. The time needed to recharge a 100% depleted battery is (5) five hours. The unit has an internal fan that ensures an even temperature throughout the cold chamber.

All units have internal and external lids, and are fully transportable with locking casters.

Displayed on the cart is: the UPS charge indicator. Displayed on the front of the refrigeration units are: the UPS battery-status, the NIST traceable data logger, and the 110V input receptacle. (The input receptacle powers both the refrigeration unit and the on-board UPS.)

VOLUME: 15 Liters

INSIDE DIMENSIONS: 11"L x 8"W x 10"H

INTERNAL BASKET: up to (12) twelve 500 ml bags of RBCs or Plasma

OUTSIDE DIMENSIONS: 31"x14"x38"

WEIGHT: 137 Lbs