



BF-15-MIL-UPSD-40AHR

15-Liter Portable Medical Refrigerator Unit with Battery Pack

GENERAL SPECIFICATIONS

VOLUME	15 Liters
Temperature Range	Refrigerator 2° to 8°C Freezer 0° to -20°C
INSIDE DIMENSIONS	11"L x 8"W x 10"H, 27.94cm x 20.32cm x 25.4cm
OUTSIDE DIMENSIONS	24" x 12" x 23", 60.96cm x 30.48cm x 58.42cm
WEIGHT	15LT Unit – 40 lbs / 18.14 kg Detachable Battery Backup System 39.5 lbs / 17.92 kg TOTAL = 79.5 lbs / 36.06 kg
Electrical Rating (AC)	100-240 VAC, 50/60Hz 0.4 Amps to 0.16 Amps
Power Consumption	60 Watts max
Compressor Capacity	18 - 150Liters Displacement 2.0cm ³
Fan	Internal fan to ensure even temperature
Datalogger (optional)	Datalogger with Temperature Probe



BATTERY

Electrical Rating (DC)	12-24 VDC ~4 Amps to 1.6Amps
Battery Life	40Ahr Battery = 26 to 32 hours
Time to Recharge	5 hours (from 100% depleted to full charge)
Battery Status	On the Cart & on the front of Refrigeration unit

BASKET (Optional)

CART (Optional)

