



## BFT-15x2

### 15-Liter Portable Medical Refrigerator Unit with Cart

The BFT-15 portable medical BioFridge with cart was designed for use within hospitals, laboratories and medical facilities. Popularly used as vaccine refrigerators and for transporting blood, it features an internal fan, sturdy handles, lockable latch and a programmable digital thermostat.

#### GENERAL SPECIFICATIONS

<b>VOLUME</b>	15 Liters
<b>Temperature Range</b>	Refrigerator 2° to 8°C Freezer 0° to -20°C
<b>INSIDE DIMENSIONS</b>	2x 11"L x 8"W x 10"H
<b>INTERNAL BASKET</b>	2x up to (12) twelve 500 ml bags of RBCs or Plasma
<b>OUTSIDE DIMENSIONS</b>	31"x14"x38"
<b>WEIGHT</b>	137 Lbs
<b>Electrical Rating (AC)</b>	100-240 VAC, 50/60Hz 0.4 Amps to 0.16 Amps
<b>Power Consumption</b>	60 Watts max
<b>Compressor Capacity</b>	18 - 150Liters Displacement 2.0cm <sup>3</sup>
<b>Fan</b>	Internal fan to ensure even temperature
<b>Datalogger (optional)</b>	Datalogger with Temperature Probe

#### BATTERY (located inside the cart)

<b>Electrical Rating (DC)</b>	12-24 VDC ~4 Amps to 1.6Amps
<b>Battery Life</b>	40Ahr Battery = 26 to 32 hours
<b>Time to Recharge</b>	5 hours (from 100% depleted to full charge)
<b>Battery Status</b>	On the Cart & on the front of Refrigeration unit

#### BASKET (Optional)

41"



34"

